100

140

Instructions for Use

Product Description:
JUC spray dressing is a water-soluble spray dressing of organosilicon quaternary ammonium salt. Organosilicon quaternary ammonium salt has the size of over 1nm, so it can't penetrate through cell membranes of microorganisms. By spraying on the body surface and forming positively-charged reticulate film, it reduces microbial load and prevents microbial penetration.

Indications:
It is indicated for protection of wounds (surgical wounds, first degree burn wounds , superficial abrasions and traumatic wounds)
It is indicated for prevention of local infections: bacterial infections (such as balanoposthitis, impetigo, etc.), Fungal infections (such as vulvoviginal candidiasis (VVC), tinea pedis, etc.), and viral infections (such as condyloma acuminatum, herpes, etc.).
It is indicated for devices such as catheters (urine catheters, ventilator catheters) for prevention of hospital-acquired infections (HAI).

Usage Method:
Clean the wound with normal saline.
Take out the bottle from the box and open the cap, press the spray head 2-7 times till the misty liquid sprays out.
Hold the bottle about 7 inches from the wound, press the spray head several times and cover the wound evenly with the liquid.
Repeat the above procedure and spray the wound evenly twice if necessary.

Dosage:
Use every 8 hours is advised.
15 days may be needed for wounds (such as surgical wounds, first degree burn wounds, superficial abrasions and traumatic wounds);
15 days may be needed for bacterial infections (such as balanoposthitis and impetigo); over 30 days may be needed for fungal infections (such as vulvovaginal candidiasis (VVC) and tinea pedis); several months may be needed for viral infections (such as condyloma acuminatum and herpes).
Long term usage for prevention of infection.

Precautions:
When pressing the spray head, the finger should press it till the end naturally for once, to avoid influencing the pulverization and dosage.

Drug Interactions:
No known drug interaction of clinical significance.

Contraindication:
JUC shall not be applied on patients who are extremely allergic to organosilicon quaternary ammonium salt.

Storage:
Storage under room temperature， shelf life 3 years.

Manufacturer:
NMS Technologies Co., Ltd.
8 Qiaobei Road, Shiqiao, Pukou District, Nanjing, 211804, China
www.juc360.com
Tel: +86 (25) 83400310
Fax: +86 (25) 83163882

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NMS Technologies Co.,Ltd.
AMT Co.,Ltd.
XI'AN GLOBAL PRINTING CO.,LTD.